



Phone : 617-800-3785 617-800-6495
Fax : 888-818-2158
Email: info@cwbiosciences.com



Ver.A/1

Disposable swab

Product name : Disposable swab

Product introduction

Disposable swab is used to collect sample of cervix, nasopharynx and oropharynx, etc. It consists of plastic rod and flocking swab head. Strict technological and process control conditions make the product free of DNA enzyme, RNA enzyme, cell inhibitor and endotoxin.

Product characteristics

The paper and plastic packaging that is convenient for sterilization is adopted.

Area of application

For clinical sampling.

Attentions

1. This product is disposable.
2. It is forbidden to use expired products or damaged inner packaging.

Instructions

1. Confirm that the product is within the validity period and the package is in good condition, tear off the package and take out the swab.
2. Sample with a swab at the desired location.

Product maintenance and transportation

1. It should be stored in a dry, ventilated room without corrosive gas.
2. General logistics transportation.

Production date

See outer packing for production date and expiration date.

Manufacturer

Jiangsu CoWin Biotech Co., Ltd.

Add: No.58 G52 Building, East of Lujia Road, West Side of Tai Road, China medical city, Taizhou, Jiangsu, 225300

Phone: 86-523-86201352

Fax: 86-523-86816890

Web: <https://www.cwbiotech.com>

E-mail: service@cwbiotech.com

EC-Representative

SUNGO Europe B.V.

VAT: NL857821659B01

Add: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Contact: SUNGO Secretary

Tel /Fax: +31 (0) 2021 11106

E-mail: ec.rep@sungogroup.com

Version: A/1

Effective Date: 2020-04-26